CUDA[®]

SURGICAL

FIBEROPTIC CABLES INSTRUCTIONS FOR USE

ST Technologies
6018 Bowdendale Avenue
Jacksonville, FL 32216

Customer Service: 904 737-7611
Toll Free 877-677-2832

RMS UK, Ltd.
28 Trinity Road
Nailsea, Somerset BS48 4NU
United Kingdom
TEL: 01275 858891

PLEASE READ BEFORE USE

Failure to follow these instructions could render the device unusable and void any warranty.

Device Description

CUDA[®] SURGICAL fiberoptic cables are designed to deliver maximum light when coupled to a medical grade fiberoptic light source. CUDA[®] SURGICAL fiberoptic cables can be used with quartz halogen, metal halide, or Xenon light sources. They are compatible with virtually all endoscopes, medical instruments, and microscopes.

Intended Use

This device is designed to illuminate a surgical site by relaying light from a fiberoptic light source onto the desired site. CUDA[®] SURGICAL cables are medical grade, high transmittance, peak efficiency.

NOTE: The style of cable suitable for a given lightsource can be determined from the table below.

| | Standard Cable | Clad Rod Cable | Lensed Cable | Fused Cable |
|-----------------------|----------------|----------------|--------------|-------------|
| Cuda 250W Halogen | X | | | |
| Other 250W Halogen | | X | | |
| Focused Ceramic Xenon | | X | | X |
| Other Xenon ≤ 300W | X | | X | |
| Xenon >300W | | | | X |
| Metal Halide ≤ 75W | X | | | |

- For use with Xenon lightsources, be sure lightsource is filtering at least 90% of infrared radiation (IR) to prevent high heat.

Contraindications

This device has no contraindications of which we have knowledge.

**Warnings and Precautions**

- CUDA[®] SURGICAL light cables are provided non-sterile and must be sterilized before use. See instructions for cleaning and sterilization.
- Always inspect cables for any evidence of damage prior to use. Pay particular attention to optical surfaces looking for scratches or dings.
- Use caution to treat CUDA[®] SURGICAL cables as you would any fine optical device.
- The user of this product should be thoroughly familiar and trained in use and care of the product.
- Match the cable size to the endoscope. Using a larger cable (fiber bundle size) than the endoscope requires will result in overheating the endoscope's proximal end.

- FIRE HAZARD: Never drape or cover the end of any fiberoptic cable with anything flammable.
- Light loss and warmer than normal instrument temperatures may be encountered if using a light source fiber optic cable with a larger aperture or bundle diameter than the instruments receiving aperture or bundle diameter.

Adverse Events

We have no knowledge of any adverse events with this product.

Cleaning of headband, headlight and headlight cable.

Exterior surface of headband and module can be wiped clean with 70% Isopropyl alcohol or mild soap. Headlight cable is not to be steam autoclaved. Autoclaving will cause damage to the non-sealed Proximal end of the cable which has a lens on it. The headlight cable can be wiped in same process as the headband and module.

Cleaning and Sterilization of Fiberoptic Light Guides

Cleaning

Fiberoptic light guides for endoscopes, microscopes, and surgical instruments are high quality optical devices. They require similar care to that taken for any precision optical component. After each use, the light guide should be washed and cleaned of all debris. A soft brush and mild detergent should be used to clean the device. Pay particular attention to any crevices or seams. Avoid any harsh materials or detergents that can scratch or in any way damage the optical surfaces on each end of the light guide. **ST Technologies** recommends the following Steam Sterilization (Autoclave) guidelines:

Steam Sterilization (Autoclave)

The cycle selected is dependent on equipment and hospital protocol. General guidelines are:

| METHOD | CYCLE | MINIMUM TEMPERATURE | EXPOSURE TIME |
|-----------------|-----------|-------------------------|---------------|
| Steam (wrapped) | Prevacuum | 270 deg. F (132 deg. C) | 4 minutes |
| Steam (wrapped) | Gravity | 270 deg. F (132 deg. C) | 10 minutes |
| Steam (wrapped) | Gravity | 250 deg. F (121 deg. C) | 45 minutes |

Chemical Disinfection

It is recommended that the guidelines below published by Advanced Sterilization, a Johnson & Johnson company be followed.

| Product | High Level Disinfection Sterilization | Sterilization |
|-------------------------------------|---|---------------|
| Cidex Activated Dialdehyde Solution | 45 minutes @ 25 degrees C 77 degrees F | 10 Hours |
| Cidex Plus 28 Day Solution | 20 minutes @ 25 degrees C 77 degrees F | NA |

Steris Sterilization

Follow Steris Corporation Guidelines

Sterrad Sterilization

Follow Sterrad (a J&J Company) Guidelines

Limited Warranty

Your fiberoptic cable has a one (1) year warranty from the date of shipment on workmanship and all defects of material except for broken fiber. Should your product prove to have such defects within one (1) year of shipment, **ST Technologies** will repair or replace the product or component part without charge.

Should your fiberoptic cable need servicing under this warranty, please contact your distributor or your customer support specialist for return authorization documentation. You should carefully pack the product in a sturdy carton, including a note describing the defects, your name, your company name, telephone number and a return address. Warranty does not cover equipment subject to misuse, accidental damage, and normal wear and tear. This warranty gives you specific legal rights, and you may also have other rights that vary from state to state.

Post warranty repairs

Please contact your distributor or your customer support specialist for return authorization documentation.

Chart of medical device symbols used

| | |
|---|--|
|  | Manufacturer |
|  | Date of manufacture (YYYYMM or YYYY) |
|  | “Authorized Representative in the European Community.” |
|  | Caution, consult accompanying documents |
|  | Caution: Federal law restricts this device to sale by or on the order of a (licensed healthcare practitioner). |
|  | Keep product dry |
|  | Batch code / lot number |
|  | CE mark |
|  | Nonsterile |
|  | Catalog, reorder or reference number |
|  | Do not use if package is damaged |

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Siga as Instruções da Steris Corporation
Contactar o seu distribuidor: EUROCASMÉDICA, S.A. – Rua do Cunha, 45 – 4200-251 Porto