CUDA[®]

SURGICAL

LED-8500II Headlight System Operation Manual

ST Technologies[®]
6018 Bowdendale Avenue
Jacksonville, FL 32216 USA

Customer Service: 904 208 2290
Toll Free 877 814 2237

RMS UK, Ltd.
28 Trinity Road
Nailsea, Somerset BS48 4NU
United Kingdom
TEL: 01275 858891

TABLE OF CONTENTS

1. INTENDED USE
2. GENERAL WARNINGS
3. ASSEMBLY
4. MAINTENANCE
5. CLEANING
6. DESKTOP SMART CHARGER AND CALIBRATOR
7. REPLACEMENT PARTS
8. USER SERVICE
9. WARRANTY AND REPAIR
10. END OF PRODUCT LIFE
11. CHART OF MEDICAL DEVICE SYMBOLS USED

1. INTENDED USE

This LED-8500II Headlight system is designed to deliver illumination from a high intensity LED for surgical site illumination.

Congratulations on the purchase of your new LED Headlight system.

The LED Headlight LED-8500II is a battery powered portable headlight with an adjustable spot and uses encrypted battery technology that requires the use of approved batteries obtained from the manufacturer or authorized distributors.

Additionally, a low battery alarm is offered in this system and will beep when power is low. Decreasing the output may extend battery life and therefore the amount of warning time before power down. This alarm is a convenience and the user should plan for backup lighting for any application using this system; see General Warnings.

2. GENERAL WARNINGS

Note:

You must charge the batteries fully before using the headlight for the first time.

The user of this product should be thoroughly familiar with the use and care of this product.

The user should carefully study the manual before making any attempt to use this product clinically.

Before each procedure, carefully check the batteries state of charge (Fuel Gauge) to assure optimal duration for the procedure.

LED High Intensity Light sources using high intensity LEDs that will produce some heat as well as light. The high output of such LEDs can cause the headlamp to feel warm to the touch during operation, this is normal.

The battery charger and the battery chargers power supply are not to be used in the operating room theater or the patient area.

Use care not to point the light directly at the eyes during operation, the brilliant light can cause severe eye discomfort.

Users are responsible for providing backup lighting systems for your application when using this device.

3. ASSEMBLING THE HEADLIGHT SYSTEM

Before continuing remove the protective lens cap.

Insert a fully charged battery into the Battery Holster and ensure it is fully seated down on the connector at the bottom of the battery cavity.

To attach the power cable from the LED Module to the Battery Holster, insert the push-pull quick connector of the LED Module into the receptacle on the lead from the Battery Holster until it clicks in place as shown here. This connection will not come off when positively engaged. You can rotate the connectors to find the keyway and then press together and not damage the pins.

The On/Off switch is also the LED intensity control. Simply rotate the knob clockwise or counter-clockwise to produce the desired light output.

The headlight has an iris to allow adjustment of the beam spot diameter. Rotate the fluted iris adjustment ring on the headlight clockwise or counterclockwise for the appropriate spot diameter.

The battery holster is worn on a belt; on the rear side of the holster is a metal belt clip that slips over and secures to a belt. The power cable secures into the cable clips on the side and rear of the headband.

To remove the LED Module from the Battery Holster, press the small raised button on the quick disconnect device located between the LED Module and the Battery Holster and pull straight out without twisting.

Next, place the headband on your head; you can make size adjustments by turning the knob on the rear of the headband until it feels comfortable. For more comfort you can adjust the cross band.

Once the headlight module and cable are connected, the unit is ready for use.

The system comes with gown clips to relieve the weight of the cable from the headband. One gown clip is attached below the shoulder.

Once the battery is depleted into the last shown “block” on the power gauge, the unit may beep once to indicate low battery levels. The unit may beep twice to indicate further reductions in power before shut-down.

4. MAINTENANCE

Storing the headlight system in a safe place will prolong the life of the unit. The following guidelines will help in sustaining the headlight system useful life:

- Keep the optical faces of the LED Headlight from touching hard surfaces, which may cause scratches on the surface. Scratches on the surface will diminish the light output.

5. CLEANING

The headlight module can be wiped down with alcohol. The lens should be cleaned only with lens tissue, available in any camera store; follow the directions on the package. Follow all applicable bloodborne pathogen procedures as required by OSHA and/or your hospital, when cleaning and disinfecting the product.

The headband can be cleaned with an alcohol wipe down or washed with a lukewarm solution of water and mild detergent. Rinse thoroughly.

Do not autoclave the LED Headlight module, headband, cable, battery or battery holster.

6. DESKTOP SMART CHARGER AND CALIBRATOR

What charger components are in the carrying bag?

1. One CH5000 desktop charger/calibrator
2. One 24V 2.5A DC power supply with universal mains input.
3. One of these main power cables.
 - North American chargers (CH5000A) are packed with a US 3-pin mains cord
 - European chargers (CH5000E) are packed with a European 2-pin mains cord
 - UK chargers (CH5000U) are packed with a UK 3-pin mains cord

Safety

1. Do not expose the charger or power supply to water or liquids; this is not a sealed case.
2. Do not open the charger or power supply case, no user serviceable parts are inside.
3. Do not cover the fan exhaust or obstruct the airflow, this will cause overheating.
4. Use only the manufacturer’s power supply and observe terminal polarity.


5. Place the charger in a cool spot, away from external heat sources

During recalibration the charger may become warm.

Using your Charger

Place the charger on a flat, level surface away from sources of heat and moisture. Plug the DC connector from the power supply into the back of the charger and connect the power supply to the mains AC supply using the cable supplied.

Charging

Place the battery into the battery bay ensuring that the 5-way connector is fully seated. The LEDs in the status window will provide status information and the charger will automatically begin charging.

LED Indication:

The status of the battery is indicated by the LEDs visible in the status window:

* Green flashing:	Battery charging
● Green solid:	Battery fully charged
* Blue flashing:	Battery in calibration mode
● Blue solid:	Battery fuel gauge calibrated
* Red flashing:	Battery fuel gauge in need of recalibration
● Red solid:	Error

Recharge Time

The ND2054 battery requires up to 3.5 hours to recharge a fully depleted battery.

Recalibration Time

Recalibration consists of a calibration charge, followed by a calibration discharge. Finally the battery is given a regular charge. A calibration cycle will be faster if the battery is fully charged to begin with. Recalibration time is governed by the battery voltage and capacity. Larger batteries will take longer to recalibrate. Calibration is initiated each time

the button is pressed, so it is not recommended to press the recalibration button part way through the recalibration cycle.

What are the SMBus and the SBS?

The Smart Battery System defines the parameters that are stored by an SBS-compliant battery. These parameters include full battery status and fuel gauging information. The System Management Bus is the language by which these parameters are communicated between the battery, the charger and the host device.

How does the charger know what charge to deliver?

The CH5000 charger is capable of sensing and delivering an appropriate charge. Upon inserting the battery into the charger, the battery communicates to the charger over the SMBus telling the charger what type of cell chemistry it is and what type of charge regime it needs. The CH5000 then configures its output to provide the charge regime requested by the battery. If no SMBus communications are issued from the battery, the charger interrogates the thermistor/resistor I.D. pin on the battery terminal and delivers an appropriate charge.

FUEL GAUGE AND RECALIBRATION:

Fuel Gauge Operation

You may receive one of two different batteries with your LED headlight system. The main difference between these two battery packs is the location and the function of the battery gauge.

On one battery type (shown below), the LCD fuel gauge employs a 'five' segment LCD. Each segment represents approximately 20% of the batteries remaining capacity. Each segment of a fully charged battery represents approximately 30 minutes of run time. For this battery, the LCD Fuel Gauge is located at the same end of the battery as the pull tab and requires no actions to view it.

On the other battery type (example shown below), the LED fuel gauge employs a 'four' segment LED. Each segment represents approximately 25% of the batteries remaining capacity. Each segment of a fully charged battery represents approximately 37.5

minutes of run time. For this battery, the LED fuel gauge is located on the front of the battery and the user is required to press the small circular button next to the gauge to determine remaining battery life, which is represented by green LEDs.

“Real-Life” applications rarely fully discharge a battery pack. Frequent partial discharges are not a problem to your smart battery, however after repeated use in this way, the accuracy of the fuel gauge may be reduced.

The smart battery has a built-in monitoring system which checks the accuracy of the fuel gauge, based on the discharge history of the battery. This history is broadcast over the SMBus and can be used by the host device to inform the user when to recalibrate the electronic fuel gauge.

Recalibration

Recalibration of the electronics is achieved by fully recharging the battery followed by a full discharge. Depending on the storage history of the battery pack, the smart battery may require calibration when brand new.

If fuel gauge recalibration is needed, the red LED on a calibrating charger will flash upon insertion of the battery. This provides feedback on the accuracy of the fuel gauge and avoids unnecessary calibration cycles.

The user can either calibrate the fuel gauge and charge the battery, or just charge the battery.

Calibration takes longer than charging and it may not be convenient to go through the calibration cycle.

To recalibrate the fuel gauge, press the button on the front of the charger.

The charger will automatically begin to charge the battery if the button is not pressed.

The blue LED will flash to indicate that the battery is undergoing the recalibration cycle.

During calibration the discharge resistors will be cooled by the fan. Removing the battery, or pressing the calibration button again will re-start the process from the beginning.

At the end of this procedure the blue LED will stay constant indicating a fully calibrated fuel gauge. Warm environments can cause calibration failure - keep the charger away from direct sunlight or heat sources.

7. REPLACEMENT PARTS

For replacement parts, additional or optional equipment :
Please call Customer Service @ 877-677-2832

8. USER SERVICE

The LED Headlight has no user or field serviceable components. It can only be serviced at the factory by factory trained technicians.

System Function and Error LEDs are visible on the back side of the Battery holster (side with the belt clip) There is one 'red' and one 'green' LED indicator lights.

When the unit is turned on:

The lit 'Green' LED signifies the system is powered up and will stay on constantly.

The lit 'Red' LED signifies the system requires service. Turn off the headlight and remove the battery. Contact Customer Service.

9. LIMITED WARRANTY

Your LED Headlight and battery holster carry a three year warranty from the date of shipment on workmanship and all defects of material.

Your batteries and charger carry a six month warranty from the date of shipment on workmanship and all defects of material. Should your product prove to have such defects **ST TECHNOLOGIES**® will repair or replace the product or component part without charge. Please contact **ST TECHNOLOGIES**® for return authorization documentation.

You should carefully pack unit in a sturdy carton and ship it to the factory. Please include a note describing the defects, your name, telephone number and a return address. Warranty does not cover equipment subject to misuse, accidental damage, normal wear and tear or if transferred to a new owner without authorization from **ST TECHNOLOGIES**®. This warranty gives you specific legal rights and you may also have other rights that vary from state to state.

POST WARRANTY REPAIRS: You may return your product(s) for repair, shipping prepaid to the factory. Your product will be inspected and an estimate of repair charges will be submitted to you for approval. Payment must be received before repairs are completed.

- Telephone: toll free 877 814 2237
- FAX number: 904 733-0012
- Customer Service: 904 208 2290

10. END OF PRODUCT LIFE

In accordance with the European Waste from Electrical and Electronic Equipment (WEEE) directive, we encourage our customers to recycle this product whenever possible. Disposal of this unit must be performed in accordance with the applicable local environmental regulations.

Each battery carries the Rechargeable Battery Recycling Corporation (RBRC) symbol and US toll-free number that can be used to locate local facilities accepting spent batteries carrying this symbol.

Recycling locations can be found at www.rbrc.com

Other recycling locations may be found at www.eiae.org

11. CHART OF MEDICAL DEVICE SYMBOLS USED

Chart of medical device symbols used

	CE mark
	Manufacturer
	“Authorized Representative in the European Community.”
	Caution, consult accompanying documents
	Do not dispose of in ordinary municipal waste, Recycle.

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