CUDA[®]

SURGICAL

LED ROVER

Instructions for Use

ST Technologies
6018 Bowdendale Avenue
Jacksonville, FL 32216
USA

Customer Service: 904 208 2291
Toll Free 877 814 2237

RMS UK, Ltd.
28 Trinity Road
Nailsea, Somerset BS48 4NU
United Kingdom
TEL: 01275 858891

TABLE OF CONTENTS

1. INTRODUCTION
2. INTENDED USE/INDICATIONS FOR USE
3. WARNINGS/LABELING
4. SPECIFICATIONS
5. OPERATION
6. BATTERY INFORMATION
7. CLEANING/STERILIZATION
8. WARRANTY
9. CHART OF MEDICAL DEVICE STMBOLS

1. INTRODUCTION

The LED Rover is a Light Emitting Diode (LED) and battery unit to couple onto the light post of commercially available endoscopes having an ACMI or Wolf light post. The LED Rover is made to take the place of conventional light cable and remote light source. It couples to the endoscope in the same fashion as a conventional light cable by press fitting onto the endoscope's light post. The LED Rover is powered by two lithium 3V, CR123 batteries (or interchangeable equivalents).

The batteries power an LED via a current converter that supplies a constant current to the LED. The light is directed into the light post of the scope. The light is then carried by fiber optics internal to the scope to illuminate the target.

2. INTENDED USE / INDICATIONS FOR USE

The LED Rover is used to illuminate the site of surgery during minimally invasive surgical procedures in arthroscopy (orthopedic surgery), laparoscopy (general and gynecological surgery) and in Endoscopy (general, gastroenterological and ENT surgery). The light is transmitted from the light source through a scope.

3. WARNINGS / LABELING

3.1 WARNINGS

 **Caution** Battery powered equipment is considered unsuitable for use in the presence of flammable anesthetics.

 **Caution** Federal law restricts this device to sale by or on the order of a licensed practitioner.

 **Caution** The LED Rover is provided **non-sterile**. If used for sterile procedures, the unit must be sterilized or draped with commercially available surgical drapes.

 **Caution** The LED Rover will become warm to the touch but even prolonged contact during use will not cause injury.













 **Caution** Avoid looking directly at the LED inside the port or pointing it directly at anyone's eyes. The LED is painfully bright, can momentarily blind a person and will cause lingering ghost images.

 **Caution** Handle damaged or leaking batteries with care. Should you come in contact with the electrolyte, wash the exposed area with soap and water. Should the eye come in contact with the electrolyte, flush the eye with water for 15 minutes and seek medical attention. Do not incinerate discarded batteries. Follow local EPA requirements for proper disposal.

 **Caution** Always follow the instructions for use provided with the rechargeable batteries and charger system which are provisioned for use with your ROVER system.

3.2 LABELING

The LED Rover will be wrapped in bubble wrap and placed in the white box labeled as shown below.

<p style="text-align: center;">CUDA[®] <i>SURGICAL</i></p> <p style="text-align: center;">ROVER III LED LIGHTSOURCE, ACMI PORT</p> <p style="text-align: center;">Caution: Federal law restricts this device to sale by or on the order of a licensed practitioner.</p> <table style="width: 100%;"> <tr> <td style="width: 15%;"></td> <td style="width: 45%;">LLS-3LBA</td> <td style="width: 40%; text-align: right;">R_x</td> </tr> <tr> <td></td> <td>WOXXXXX</td> <td></td> </tr> <tr> <td></td> <td>ST Technologies 6018 Bowdendale Ave. Jacksonville, FL 32216 U.S.A.</td> <td></td> </tr> <tr> <td></td> <td>20YY</td> <td style="text-align: right;">CE</td> </tr> <tr> <td></td> <td>RMS UK LTD 28 Trinity Rd. Nailsea, Somerset BS48 4NU United Kingdom</td> <td></td> </tr> </table>		LLS-3LBA	R_x		WOXXXXX			ST Technologies 6018 Bowdendale Ave. Jacksonville, FL 32216 U.S.A.			20YY	CE		RMS UK LTD 28 Trinity Rd. Nailsea, Somerset BS48 4NU United Kingdom		<p style="text-align: center;">CUDA[®] <i>SURGICAL</i></p> <p style="text-align: center;">ROVER III LED LIGHTSOURCE, ACMI PORT RECHARGEABLE BATTERIES</p> <p style="text-align: center;">Caution: Federal law restricts this device to sale by or on the order of a licensed practitioner.</p> <table style="width: 100%;"> <tr> <td style="width: 15%;"></td> <td style="width: 45%;">LLS-3LBA-R</td> <td style="width: 40%; text-align: right;">R_x</td> </tr> <tr> <td></td> <td>WOXXXXX</td> <td></td> </tr> <tr> <td></td> <td>ST Technologies 6018 Bowdendale Ave. Jacksonville, FL 32216 U.S.A.</td> <td></td> </tr> <tr> <td></td> <td>20YY</td> <td style="text-align: right;">CE</td> </tr> <tr> <td></td> <td>RMS UK LTD 28 Trinity Rd. Nailsea, Somerset BS48 4NU United Kingdom</td> <td></td> </tr> </table>		LLS-3LBA-R	R_x		WOXXXXX			ST Technologies 6018 Bowdendale Ave. Jacksonville, FL 32216 U.S.A.			20YY	CE		RMS UK LTD 28 Trinity Rd. Nailsea, Somerset BS48 4NU United Kingdom	
	LLS-3LBA	R_x																													
	WOXXXXX																														
	ST Technologies 6018 Bowdendale Ave. Jacksonville, FL 32216 U.S.A.																														
	20YY	CE																													
	RMS UK LTD 28 Trinity Rd. Nailsea, Somerset BS48 4NU United Kingdom																														
	LLS-3LBA-R	R_x																													
	WOXXXXX																														
	ST Technologies 6018 Bowdendale Ave. Jacksonville, FL 32216 U.S.A.																														
	20YY	CE																													
	RMS UK LTD 28 Trinity Rd. Nailsea, Somerset BS48 4NU United Kingdom																														
<p style="text-align: center;">CUDA[®] <i>SURGICAL</i></p> <p style="text-align: center;">ROVER III LED LIGHTSOURCE, WOLF PORT</p> <p style="text-align: center;">Caution: Federal law restricts this device to sale by or on the order of a licensed practitioner.</p> <table style="width: 100%;"> <tr> <td style="width: 15%;"></td> <td style="width: 45%;">LLS-3LBW</td> <td style="width: 40%; text-align: right;">R_x</td> </tr> <tr> <td></td> <td>WOXXXXX</td> <td></td> </tr> <tr> <td></td> <td>ST Technologies 6018 Bowdendale Ave. Jacksonville, FL 32216 U.S.A.</td> <td></td> </tr> <tr> <td></td> <td>20YY</td> <td style="text-align: right;">CE</td> </tr> <tr> <td></td> <td>RMS UK LTD 28 Trinity Rd. Nailsea, Somerset BS48 4NU United Kingdom</td> <td></td> </tr> </table>		LLS-3LBW	R_x		WOXXXXX			ST Technologies 6018 Bowdendale Ave. Jacksonville, FL 32216 U.S.A.			20YY	CE		RMS UK LTD 28 Trinity Rd. Nailsea, Somerset BS48 4NU United Kingdom		<p style="text-align: center;">CUDA[®] <i>SURGICAL</i></p> <p style="text-align: center;">ROVER III LED LIGHTSOURCE, WOLF PORT RECHARGEABLE BATTERIES</p> <p style="text-align: center;">Caution: Federal law restricts this device to sale by or on the order of a licensed practitioner.</p> <table style="width: 100%;"> <tr> <td style="width: 15%;"></td> <td style="width: 45%;">LLS-3LBW-R</td> <td style="width: 40%; text-align: right;">R_x</td> </tr> <tr> <td></td> <td>WOXXXXX</td> <td></td> </tr> <tr> <td></td> <td>ST Technologies 6018 Bowdendale Ave. Jacksonville, FL 32216 U.S.A.</td> <td></td> </tr> <tr> <td></td> <td>20YY</td> <td style="text-align: right;">CE</td> </tr> <tr> <td></td> <td>RMS UK LTD 28 Trinity Rd. Nailsea, Somerset BS48 4NU United Kingdom</td> <td></td> </tr> </table>		LLS-3LBW-R	R_x		WOXXXXX			ST Technologies 6018 Bowdendale Ave. Jacksonville, FL 32216 U.S.A.			20YY	CE		RMS UK LTD 28 Trinity Rd. Nailsea, Somerset BS48 4NU United Kingdom	
	LLS-3LBW	R_x																													
	WOXXXXX																														
	ST Technologies 6018 Bowdendale Ave. Jacksonville, FL 32216 U.S.A.																														
	20YY	CE																													
	RMS UK LTD 28 Trinity Rd. Nailsea, Somerset BS48 4NU United Kingdom																														
	LLS-3LBW-R	R_x																													
	WOXXXXX																														
	ST Technologies 6018 Bowdendale Ave. Jacksonville, FL 32216 U.S.A.																														
	20YY	CE																													
	RMS UK LTD 28 Trinity Rd. Nailsea, Somerset BS48 4NU United Kingdom																														

The device has no known contraindications.

4. SPECIFICATIONS

Item	Specification
Lamp Type	LED
Color temperature	6500° K (nominal)
Lamp life	50,000 hrs. (nominal)
Battery	Lithium CR123, 3Volts (2X) or Lithium Ion RCR123A 3V (2X)
Current converter output	611 milliamps (nominal)
Discharge time	110 minutes (typical, actual endurance may vary)
Model	LED Rover
Dimensions	5.1 inches long X 1.25 inches diameter at LED housing
Weight	3.7 oz.
Electrical Protection	Class I, Type BF
Medical Classification	FDA Class II Europe Class I

5. OPERATION

The LED Rover is operated by a pushbutton switch inside a protective silicone rubber boot at the back end of the handle. It features a momentary ON function with partial depression of the switch, and clicks to stay ON with full depression. Depress the switch fully again to turn it off.

6. BATTERY INFORMATION

6.1 DISPOSABLE BATTERY MODEL BATTERY REPLACEMENT

Unscrew the headpiece that houses the LED. (Do not remove the switch; the batteries cannot fit through the switch end.) Remove the old batteries and insert the new ones. Observe correct polarity. The positive (+) end of both batteries must point toward the headpiece. The positive end of the battery has a protruding nipple, the negative (-) end is flat.

NOTE: *Because discharge characteristics may vary between designations or between used and unused batteries, always replace batteries in pairs with pairs of identical new batteries.*

Replace disposable 3 volt Lithium CR123 primary batteries only with similar size and voltage batteries.

Storage: Store in a dry area at 40° F – 85° F.

6.2 RECHARGEABLE BATTERY MODEL

The rechargeable model lightsource differs from the standard only by the use of rechargeable batteries.

The batteries supplied may be of various voltages which work in the ROVER system.

All batteries will be RCR123R-type cells.

Always use batteries of the same chemistries and voltages together in the ROVER system.

Battery life in the ROVER system will vary depending upon the battery voltages used.

Always replace batteries in pairs which have each been fully charged.

Always replace batteries in pairs when discarding units that are spent beyond their useful life. Follow all local and government policies for battery disposal.

Consult the instructions for use with your rechargeable battery system for proper operation and warning statements.

Never use chargers that are not designed to work with the battery chemistry and voltages of cells supplied with your ROVER system.

Battery recharge times will vary based upon the battery voltages used and the charger system supplied.

7. CLEANING / STERILIZATION

7.1 CLEANING

The LED Rover is a high quality optical device. Care for it appropriately.

A soft brush with mild detergent may be used to clean the device. The lens of the LED inside the optic port may be wiped with a swab. Avoid harsh or abrasive material that could scratch or in any way damage the optical surface.

7.2 STERILIZATION

The LED Rover is provided non-sterile. The device should be cleaned and disinfected or sterilized prior to use, according to your institution's requirements.

Caution ***Do not reprocess the LED Rover using steam sterilization or autoclave. Use of these processes will result in damage to the instrument, and void its warranty.***

Process per your institution's validated procedures and process parameters.

Compatible Sterilization Method: STERRAD™ in STERRAD sterilizer. Double wrap with sterilization wrap. The full cycle of the Sterrad sterilizer comprises two identical exposure phases (6 minute injection, 2 minute diffusion and 2 minute plasma) preceded by a 10-minute pre-exposure plasma phase.

Follow the instructions provided with the STERRAD Sterilization System for sterilization processing via this equipment and method.

NOTE: The LED Rover has a black anodized exterior finish which is affected by Sterrad sterilization. The black anodized finish gradually bleaches; however, this is irrelevant to the functionality of the unit.

7.3 HIGH LEVEL DISINFECTION (IMMERSION)


The LED Rover may be immersed in high level disinfectant solution no deeper than six (6) inches (15cm) for as long as one (1) hour. Process according to the disinfectant manufacturer's instructions or your institution's established protocol. We make no claims other than the Rover's immersibility as stated.

8. LIMITED WARRANTY

Your fiberoptic products(s) carry a 1-year warranty from the date of shipment on workmanship and all defects of material except batteries. Should your product prove to have such defects within one year from the date of shipment, **ST TECHNOLOGIES** will repair or replace the product or component part without charge. Should your fiberoptic product(s) need servicing under this warranty, please contact **ST TECHNOLOGIES** for return authorization documentation. You should carefully pack product in a sturdy carton and ship it to the factory. Please include a note describing the defects, your name, telephone number and a return address. Warranty does not cover equipment subject to misuse, accidental damage, and normal wear and tear. This warranty gives you specific legal rights and you may also have other rights that vary from state to state.

- In the U.S call: 877 814 2237 (toll free)
- Customer Service 904 208 2291

9. CHART OF MEDICAL DEVICE SYMBOLS USED

	Manufacturer
	Date of manufacture (YYYYMM or YYYY)
	“Authorized Representative” in the European Community
	Caution, consult accompanying documents
	Consult Instructions for Use
	CE mark

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
<p style="text-align: center;">CUDA[®] <i>SURGICAL</i></p> <p style="text-align: center;">ROVER III LED LIGHTSOURCE, ACMI PORT</p> <p style="text-align: center;"><small>Caution: Federal law restricts this device to sale by or on the order of a licensed practitioner.</small></p> <p>REF LLS-3LBA LOT WOXXXXX</p> <p> ST Technologies 6018 Bowdendale Ave. Jacksonville, FL 32216 U.S.A.</p> <p> 20YY</p> <p> RMS UK LTD 28 Trinity Rd. Nailsea, Somerset BS48 4NU United Kingdom</p> <p style="text-align: right;">R_x CE</p>	<p style="text-align: center;">CUDA[®] <i>SURGICAL</i></p> <p style="text-align: center;">ROVER III LED LIGHTSOURCE, ACMI PORT RECHARGEABLE BATTERIES</p> <p style="text-align: center;"><small>Caution: Federal law restricts this device to sale by or on the order of a licensed practitioner.</small></p> <p>REF LLS-3LBA-R LOT WOXXXXX</p> <p> ST Technologies 6018 Bowdendale Ave. Jacksonville, FL 32216 U.S.A.</p> <p> 20YY</p> <p> RMS UK LTD 28 Trinity Rd. Nailsea, Somerset BS48 4NU United Kingdom</p> <p style="text-align: right;">R_x CE</p>
<p style="text-align: center;">CUDA[®] <i>SURGICAL</i></p> <p style="text-align: center;">ROVER III LED LIGHTSOURCE, WOLF PORT</p> <p style="text-align: center;"><small>Caution: Federal law restricts this device to sale by or on the order of a licensed practitioner.</small></p> <p>REF LLS-3LBW LOT WOXXXXX</p> <p> ST Technologies 6018 Bowdendale Ave. Jacksonville, FL 32216 U.S.A.</p> <p> 20YY</p> <p> RMS UK LTD 28 Trinity Rd. Nailsea, Somerset BS48 4NU United Kingdom</p> <p style="text-align: right;">R_x CE</p>	<p style="text-align: center;">CUDA[®] <i>SURGICAL</i></p> <p style="text-align: center;">ROVER III LED LIGHTSOURCE, WOLF PORT RECHARGEABLE BATTERIES</p> <p style="text-align: center;"><small>Caution: Federal law restricts this device to sale by or on the order of a licensed practitioner.</small></p> <p>REF LLS-3LBW-R LOT WOXXXXX</p> <p> ST Technologies 6018 Bowdendale Ave. Jacksonville, FL 32216 U.S.A.</p> <p> 20YY</p> <p> RMS UK LTD 28 Trinity Rd. Nailsea, Somerset BS48 4NU United Kingdom</p> <p style="text-align: right;">R_x CE</p>

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<p style="text-align: center;">CUDA[®] <i>SURGICAL</i></p> <p style="text-align: center;">ROVER III LED LIGHTSOURCE, ACMI PORT Caution: Federal law restricts this device to sale by or on the order of a licensed practitioner.</p> <p>REF LLS-3LBA LOT WOXXXX</p> <p> ST Technologies 6018 Bowdendale Ave. Jacksonville, FL 32216 U.S.A.</p> <p> 20YY</p> <p>EC REP RMS UK LTD 28 Trinity Rd. Nailsea, Somerset BS48 4NU United Kingdom</p> <p style="text-align: right;">R_x CE</p>	<p style="text-align: center;">CUDA[®] <i>SURGICAL</i></p> <p style="text-align: center;">ROVER III LED LIGHTSOURCE, ACMI PORT RECHARGEABLE BATTERIES Caution: Federal law restricts this device to sale by or on the order of a licensed practitioner.</p> <p>REF LLS-3LBA-R LOT WOXXXX</p> <p> ST Technologies 6018 Bowdendale Ave. Jacksonville, FL 32216 U.S.A.</p> <p> 20YY</p> <p>EC REP RMS UK LTD 28 Trinity Rd. Nailsea, Somerset BS48 4NU United Kingdom</p> <p style="text-align: right;">R_x CE</p>
<p style="text-align: center;">CUDA[®] <i>SURGICAL</i></p> <p style="text-align: center;">ROVER III LED LIGHTSOURCE, WOLF PORT Caution: Federal law restricts this device to sale by or on the order of a licensed practitioner.</p> <p>REF LLS-3LBW LOT WOXXXX</p> <p> ST Technologies 6018 Bowdendale Ave. Jacksonville, FL 32216 U.S.A.</p> <p> 20YY</p> <p>EC REP RMS UK LTD 28 Trinity Rd. Nailsea, Somerset BS48 4NU United Kingdom</p> <p style="text-align: right;">R_x CE</p>	<p style="text-align: center;">CUDA[®] <i>SURGICAL</i></p> <p style="text-align: center;">ROVER III LED LIGHTSOURCE, WOLF PORT RECHARGEABLE BATTERIES Caution: Federal law restricts this device to sale by or on the order of a licensed practitioner.</p> <p>REF LLS-3LBW-R LOT WOXXXX</p> <p> ST Technologies 6018 Bowdendale Ave. Jacksonville, FL 32216 U.S.A.</p> <p> 20YY</p> <p>EC REP RMS UK LTD 28 Trinity Rd. Nailsea, Somerset BS48 4NU United Kingdom</p> <p style="text-align: right;">R_x CE</p>

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<p style="text-align: center;">CUDA[®] <i>SURGICAL</i></p> <p style="text-align: center;">ROVER III LED LIGHTSOURCE, ACMI PORT</p> <p style="text-align: center;"><small>Caution: Federal law restricts this device to sale by or on the order of a licensed practitioner.</small></p> <table style="width: 100%; border: none;"> <tr> <td style="border: 1px solid black; padding: 2px;">REF</td> <td style="padding: 2px;">LLS-3LBA</td> <td rowspan="2" style="font-size: 2em; vertical-align: middle; text-align: center;">Rx</td> </tr> <tr> <td style="border: 1px solid black; padding: 2px;">LOT</td> <td style="padding: 2px;">WOXXXXX</td> </tr> </table> <p style="font-size: 0.8em;">ST Technologies 6018 Bowdendale Ave. Jacksonville, FL 32216 U.S.A.</p> <p style="font-size: 0.8em;">20YY</p>	REF	LLS-3LBA	Rx	LOT	WOXXXXX	CE
REF	LLS-3LBA	Rx				
LOT	WOXXXXX					
EC REP	RMS UK LTD 28 Trinity Rd. Nailsea, Somerset BS48 4NU United Kingdom					
<p style="text-align: center;">CUDA[®] <i>SURGICAL</i></p> <p style="text-align: center;">ROVER III LED LIGHTSOURCE, ACMI PORT RECHARGEABLE BATTERIES</p> <p style="text-align: center;"><small>Caution: Federal law restricts this device to sale by or on the order of a licensed practitioner.</small></p> <table style="width: 100%; border: none;"> <tr> <td style="border: 1px solid black; padding: 2px;">REF</td> <td style="padding: 2px;">LLS-3LBA-R</td> <td rowspan="2" style="font-size: 2em; vertical-align: middle; text-align: center;">Rx</td> </tr> <tr> <td style="border: 1px solid black; padding: 2px;">LOT</td> <td style="padding: 2px;">WOXXXXX</td> </tr> </table> <p style="font-size: 0.8em;">ST Technologies 6018 Bowdendale Ave. Jacksonville, FL 32216 U.S.A.</p> <p style="font-size: 0.8em;">20YY</p>	REF	LLS-3LBA-R	Rx	LOT	WOXXXXX	CE
REF	LLS-3LBA-R	Rx				
LOT	WOXXXXX					
EC REP	RMS UK LTD 28 Trinity Rd. Nailsea, Somerset BS48 4NU United Kingdom					
<p style="text-align: center;">CUDA[®] <i>SURGICAL</i></p> <p style="text-align: center;">ROVER III LED LIGHTSOURCE, WOLF PORT</p> <p style="text-align: center;"><small>Caution: Federal law restricts this device to sale by or on the order of a licensed practitioner.</small></p> <table style="width: 100%; border: none;"> <tr> <td style="border: 1px solid black; padding: 2px;">REF</td> <td style="padding: 2px;">LLS-3LBW</td> <td rowspan="2" style="font-size: 2em; vertical-align: middle; text-align: center;">Rx</td> </tr> <tr> <td style="border: 1px solid black; padding: 2px;">LOT</td> <td style="padding: 2px;">WOXXXXX</td> </tr> </table> <p style="font-size: 0.8em;">ST Technologies 6018 Bowdendale Ave. Jacksonville, FL 32216 U.S.A.</p> <p style="font-size: 0.8em;">20YY</p>	REF	LLS-3LBW	Rx	LOT	WOXXXXX	CE
REF	LLS-3LBW	Rx				
LOT	WOXXXXX					
EC REP	RMS UK LTD 28 Trinity Rd. Nailsea, Somerset BS48 4NU United Kingdom					
<p style="text-align: center;">CUDA[®] <i>SURGICAL</i></p> <p style="text-align: center;">ROVER III LED LIGHTSOURCE, WOLF PORT RECHARGEABLE BATTERIES</p> <p style="text-align: center;"><small>Caution: Federal law restricts this device to sale by or on the order of a licensed practitioner.</small></p> <table style="width: 100%; border: none;"> <tr> <td style="border: 1px solid black; padding: 2px;">REF</td> <td style="padding: 2px;">LLS-3LBW-R</td> <td rowspan="2" style="font-size: 2em; vertical-align: middle; text-align: center;">Rx</td> </tr> <tr> <td style="border: 1px solid black; padding: 2px;">LOT</td> <td style="padding: 2px;">WOXXXXX</td> </tr> </table> <p style="font-size: 0.8em;">ST Technologies 6018 Bowdendale Ave. Jacksonville, FL 32216 U.S.A.</p> <p style="font-size: 0.8em;">20YY</p>	REF	LLS-3LBW-R	Rx	LOT	WOXXXXX	CE
REF	LLS-3LBW-R	Rx				
LOT	WOXXXXX					
EC REP	RMS UK LTD 28 Trinity Rd. Nailsea, Somerset BS48 4NU United Kingdom					

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