# CUDA<sup>®</sup>

## SURGICAL

# Elite 1000

(Patent Pending)

## Xenon Illuminator with Manual Dimming Control Operator Manual

ST Technologies<sup>®</sup>  
6018 Bowdendale Avenue  
Jacksonville, FL 32216 USA

Customer Service: 904 737 7611  
Toll Free 877 677 2832

RMS UK, Ltd.  
28 Trinity Road  
Nailsea, Somerset BS48 4NU  
United Kingdom  
TEL: 01275 858891

## TABLE OF CONTENTS

1. INTRODUCTION
2. WARNINGS
3. SPECIFICATIONS
4. OPERATING ELEMENTS, SYMBOLS AND FUNCTIONS
  - 4.1 FRONT PANEL
  - 4.2 REAR PANEL
5. INSTALLATION
  - 5.1 SETTING UP LIGHTSOURCE
  - 5.2 CONNECTING LIGHT CABLE
6. OPERATION
  - 6.1 POWERING UP LIGHTSOURCE
  - 6.2 LIGHT BRIGHTNESS CONTROL
  - 6.3 LAMP LIFE HOURS DISPLAY
7. CLEANING and DISINFECTION
8. MAINTENANCE, SERVICING & REPAIR, WARRANTY
  - 8.1 LAMP REPLACEMENT
  - 8.2 WARRANTY
9. END OF PRODUCT LIFE
10. TROUBLESHOOTING
11. CHART OF MEDICAL DEVICE SYMBOLS USED
12. CHART OF ELECTRICAL SYMBOLS USED

## 1. INTRODUCTION

Congratulations on the purchase of your new Elite 1000 *Xenon Illuminator!*

This user-friendly Xenon Illuminator is a high efficiency light source utilizing state-of-the-art illumination technology. It offers a variety of features such as:

- 5600 K daylight brightness for perfect color definition
- Quiet operation
- Auto shutter blocks light when cable removed
- Turret which adapts to various types of light guides
- Mechanical dimming
- Easy lamp replacement with precision mechanical alignment
- Hour Meter (resettable)

These operator instructions for Use Manual will help you to install the device and optimally integrate it with other components of your system. It will instruct you on how to operate the Xenon Illuminator and how to keep it clean. It will give you maintenance and service guidelines as well as recommendations for best performance results.

### **INTENDED USE**

The intended use of this device is to provide light for fiberoptic cables and instruments – providing light for instrumentation via fiberoptic cables for use in surgical fields.

## 2. WARNINGS/ CAUTION

 **WARNING** The illuminator uses a highly concentrated light source (luminous power per area) and this high energy density is retained through connected lightguides and instruments. The output of a connected instrument left in close proximity or contact with tissue or flammable materials presents a risk of injury or fire. Qualified personnel must determine a safe working distance and intensity setting for each application. The output should never be left on unattended.

 **Caution** Rx only. Federal law restricts this device to sale by or on the order of a licensed healthcare practitioner.

 **Caution** To prevent fire or electric shock, do not open or expose the lightsource unit to rain or moisture. Refer all servicing to qualified personnel only.

 **Caution** Not suitable for use in presence of flammable anesthetic mixture with air or with oxygen or nitrous oxide.

 **Caution** To prevent any potential electro-magnetic interference, do not use any kind of cellular phone near the lightsource.

 **Caution** This product should be used only with type BF endoscopic instruments which have been certified according to IEC 60601-1 for medical equipment and IEC 60601-2-18 for endoscopic equipment.

This symbol indicates type BF equipment.

 **Caution** User must not alter this device in any fashion. Doing so voids all warranties and statements of suitability for any purpose.

 **Caution** All devices connecting to the Illuminator must be classified as medical equipment. Additional information processing equipment connected to the Illuminator, a Medical System and the operator must determine that all equipment complies with the appropriate end-product standards (such as IEC 60950 or IEC 60065 and the Standard for Medical System, IEC 60601-1-1).

 **Caution** Always set the intensity control to the minimum level and insert the fiberoptic cable into the unit before turning **on** the power. When light is not required at the surgical site, the intensity control should be set to the fully dimmed position. If it becomes necessary to remove the fiberoptic cable without turning the unit **off**, turn the intensity control to the fully dimmed position.

 **Caution** **The fiberoptic cable must be electrically NON-CONDUCTIVE.** It should not have conductive shielding or any other conductive connection between the patient and equipment. Such connection will impair safety of the equipment. It must be rinsed free of soaking/disinfectant solution and dried before plugging into the lightsource receptacle. Ensure the optical surface is clean before connecting the cable to the lightsource.

 **Caution** When using high intensity lightsources at full output, the recommended distance from the headlight module to the patient is not less than 12 inches (30.5 cm). If using less than 12 inches from the patient, the light intensity must be reduced.

 **Caution** User is responsible for determining if interruption of light output will create an unacceptable risk. If this determination is made arrangements (such as a standby lightsource) should be made to reduce the risk.

 **Caution** Use of this illuminator with accessories or attachments requires that the end-user follow all accessory or attachment instructions which could affect illuminator setup, usage and/or settings.

 **Caution** This device contains stored electrical energy, follow operating instructions when changing lamp.

 **Caution** Do not remove lamp immediately after operation. Allow lamp to cool 10 minutes.

 **Caution** To prevent overheating, replace only with the same type and rating of lamp. Read instructions before replacing lamp.

 **Caution** Do not operate lightsource without lamp in place.

 **Caution** For grounding reliability, use only hospital grade marked power cord and receptacle.

 **Caution** Disconnect power supply cord before servicing to avoid electrical shock.

 **Caution** To reduce risk of electrical shock, do not remove cover. Refer servicing to qualified personal.

 **Caution** Lightsource can cause permanent eye damage if viewed directly with unprotected eye.

 **Caution** Do not reach hand inside the unit due to risk of shock hazard.

### 3. SPECIFICATIONS

Item	Specification
Lamp type	Ceramic type 300 watt Xenon
Power	300 Watt
Color temperature	5600 K
Lamp life	1000 hours (typical)
Lamp replacement	Cartridge replacement from the side of the unit
Light guide adapter	Turret type to accommodate standard Storz, ACMI, Wolf, Olympus cable endtips
Brightness control	Mechanical progressive shutter control from the front panel
Input voltage	100-120 V AC, 50/60 Hz 220-240 V AC, 50 Hz
Power consumption	450 watt max
Regulatory Approvals	UL60601, EN 60601-1, EN 60601-1-2, CAN/CSA, C22.2, No 601.1-M90
Equipment Class	BF-type
Mode of Operation	Continuous operation
Water Resistant	Not Protected Equipment, IPX0
Operating Environment Temperature Relative Humidity Air Pressure	+20° to +40° C (68° to 104° F) 30 to 85% 700 to 1060 hPa
Storage Environment Temperature Relative Humidity Air Pressure	-20° to +60° C (-4° to 140° F) 0 to 95% 700 to 1060 hPa
Dimensions, mm (inch)	338 (13.3)Width x 155 (6.1)Height x 457 (18)Depth front-to-back
Weight, kg (lb)	10kg (22lb)

## 4. OPERATING ELEMENTS, SYMBOLS AND FUNCTIONS

### 4.1 FRONT PANEL


Figure 1. Illuminator Front Panel

No.	Name	Function
1.	Lamp switch	Turns the illuminator on and off; illuminates blue when turned on
2.	Intensity Control	Mechanically controls the light output
3.	Turret	Auto shutter turret accepts end tip of fiber optic cable. Light is shuttered off when no cable is present.

## 4.2 REAR PANEL


Figure 2. Lightsource Rear Panel

No.	Name	Function
1.	Fan	Fan 12V DC, airflow for cooling of unit
2.	Circuit breaker	Over current protection
3.	AC mains power switch and input module	Accepts AC Power cord and switches on AC power
4.	Hour Meter	Allows tracking of total running time on unit and provides a control on individual lamp hours
5.	Grounding Connector	For electric potential equalization
6.	Product label	Product Label containing: Model Number, Serial Number, Part Number, Electrical Ratings, Manufacturer Name and Date, Regulatory Marks, UDI Information and FDA "Rx Only" Symbol.

## 5. INSTALLATION

### 5.1 SETTING UP THE LIGHTSOURCE

Place the illuminator on a stable surface (cart, counter, stand, etc.).

**NOTE** Avoid places where the illuminator may be splashed with liquid. Absolutely **DO NOT** use in any environment with explosive or flammable gases.

**This illuminator pulls cooling air into its cabinet through a filtered grille in the bottom panel underneath the front portion and exhausts warm air from the fan in the rear.**

**DO NOT block the space in front or below the front of the illuminator or impede exhaust air flow behind or above the cabinet.**

**A clear space minimum of 5"/12.7cm behind and above the cabinet is required.**

**The illuminator should not be placed where its exhaust will influence other devices, nor where the heated air from other devices will influence the illuminator.**

Ensure the rear power switch is in OFF position.

Ensure that the lamp is properly positioned for use. Open the lamp compartment door. Move the lamp positioning lever (found to the right of the lamp) clockwise, from horizontal to vertical. As the positioning cones disengage the lamp should drop about 3mm. Move the lamp positioning lever counterclockwise to reengage the positioning cones and observe that the lamp rises about 3mm as it returns to operating position. Close the lamp compartment door.

Connect AC power cord to the power inlet located on the rear panel of the illuminator.

Ensure the power cord is completely seated into the power inlet and engage the cord retainer if so equipped.

**Use only power cords provided with the illuminator.**

Plug the AC power cord into a wall outlet using the power cord plug supplied with the unit.

**To prevent electric shock, connect power cords of peripheral equipment through medical isolation transformers.**

**NOTE**

When using medical isolation transformer, be sure to check the transformer power ratings.

Ensure that the power cord is connected to the main power with three prong plug (USA use UL2601-1 rated isolation transformers and/or power strips only).

**CAUTION**

Nearby radio or TV units may subject this equipment to radio interference. To avoid adverse electromagnetic effects, DO NOT operate this equipment near RF energy equipment.

## 5.2 CONNECTING THE LIGHT CABLE

Connect the fiberoptic cable to endoscope or other instrument and then plug the cable illuminator end tip into the appropriate port of the turret on the front panel.

## 6. OPERATION

### 6.1 POWERING UP THE LIGHTSOURCE

To operate the illuminator: First attach the illuminator to the instrument or headlamp.

Turn on the rear power switch. Turn on the lamp switch. The fan will start. The indicator on lamp switch will light.

### 6.2 LIGHT BRIGHTNESS DISPLAY AND CONTROL

The light brightness adjustment enables the user to obtain a brighter or dimmer illumination of the object of observation.

Adjust the light intensity by sliding the intensity control knob from left (off) to right (brightest) as shown by the ramp graphics on the dimmer control.

**NOTE**

The brightness control operates a progressive mechanical shutter. The lamp is running at full power if the power switch is on no matter where the brightness control is set.

### 6.3 LAMP LIFE HOURS DISPLAY

The hour meter shows the total elapsed running time of the unit.

## 7. CLEANING and DISINFECTION

**NOTE**

***Always disconnect the power cord before cleaning the system and while unit is drying if wet-wiped.***

Disinfection agents should be commercially available cleansers commonly used for disinfection of electronic equipment cabinets in hospitals, such as ethyl or isopropyl

alcohols, disinfecting sprays containing quaternary ammonium compounds or hydrogen peroxide.

Do not use strongly caustic or acidic cleansers such as “Clorox” hypochlorite bleach, ammonia, muriatic acid or similar products. Do not use acetone, methyl ethyl ketone or halogenated / chlorinated hydrocarbon solvents or cleaners containing any of these restricted compounds.

Apply cleaning agents by light spray or dampened towels. Do not pour liquids over the cabinet.

Do not allow frank liquids to enter the cabinet seams or ventilation openings.

## 8. MAINTENANCE, SERVICING AND REPAIR, WARRANTY

Performance of preventive maintenance is not essential. Regular maintenance can, however, contribute to identifying potential problems before they become serious; thus, enhancing the instrument’s reliability and extending its useful operating life. Maintenance services can be obtained from your local representative or from the manufacturer.

Defective items or equipment are to be serviced and repaired exclusively by persons authorized by the manufacturer. All repair work shall employ original manufacturer’s parts only.

### 8.1 LAMP REPLACEMENT

<b>WARNING</b>	
<b><i>CERAMIC lamps are at high internal pressure when cold and at operating temperature; therefore, CERAMIC lamps may unexpectedly rupture resulting in discharge of hot fragments of quartz and/or glass and metal.</i></b>	
<b><i>Only handle lamps with protective covers in place. Do not handle lamps without their protective covers unless government approved safety glasses, facemasks (with neck protector, chest protector and gauntlets) are worn.</i></b>	

 <b>CAUTION</b>	<b><i>Always disconnect power cord and turn main switch off before lamp replacement.</i></b>
---	--

 <b>CAUTION</b>	<b><i>Take utmost care not to hit gas-seal edge of lamp; it may cause lamp failure.</i></b>
---	---

 <b>CAUTION</b>	<b><i>Do not touch the lamp and patient simultaneously when changing the lamp.</i></b>
---	--

 <b>CAUTION</b>	<b><i>Lamp power connection pins at the back of the lamp compartment may retain stored energy. Do not contact the lamp power connectors in the back of the lamp compartment.</i></b>
---	--

**This illuminator is designed for exclusive use with the ST TECHNOLOGIES p/n SYC0300 lamp module. Failure to use this lamp could void the warranty.**

Open the lamp compartment door.

Move the lamp positioning lever (found to the right of the lamp) clockwise, from horizontal to vertical, and pull the lamp cartridge out.

Grasp only the horizontal top and bottom grip flanges of the plastic housing of the lamp, pull straight out to remove the lamp. Insert a new one (lamp pointed to the left), pushing firmly to assure full engagement to the power connectors. Lower the lamp engagement lever counter-clockwise, from vertical to horizontal, and observe that the lamp rises about 3mm as it moves to its proper operating position.

Log the serial number and hours on the hour meter. Close the lamp compartment door. Re-connect the power cord and turn the light source on according to 6.1. Reset hour meter to zero by depressing hour meter reset button on rear panel.

<b>NOTE</b>	If the tamper seal on the lightsource is broken, the WARRANTY WILL BE VOIDED. There are no other user serviceable parts in the unit.
-------------	--

## 8.2 LIMITED WARRANTY

Your illumination product(s) carry(s) a 3-year warranty from the date of shipment on workmanship and all defects of material, excluding, fiberoptic cable and lamps. Should your product prove to have such defects within three years of the shipment, **ST Technologies**® will repair or replace the product or component part without charge. Should your illumination product(s) need servicing under this warranty, please contact **ST Technologies**® for return authorization documentation. Please carefully pack the unit in a sturdy carton and ship it to the factory. Include a note describing the defects, your name, telephone number and a return address. Warranty does not cover equipment subject to misuse, accidental damage, normal wear and tear or if transferred to a new owner without authorization from **ST Technologies**®. This warranty gives you specific legal rights and you may also have other rights that vary from state to state.

**POST WARRANTY REPAIRS:** You may return your product(s) for repair, shipping prepaid to the factory. Your product will be inspected and an estimate of repair charges will be submitted to you for approval.

- In the US call: 877 677-2832 (toll free)
- FAX number: 904 733-4832

## 9. END OF PRODUCT LIFE

In accordance with the European Waste from Electrical and Electronic Equipment (WEEE) directive, we encourage our customers to recycle this product whenever possible. Disposal of this unit must be performed in accordance with the applicable local environmental regulations.

In the US a list of recyclers in your area can be found at: <http://www.eiae.org/>.

Please contact customer service to issue a return authorization to return product to manufacturer at the end of product life.

## 10. TROUBLESHOOTING

<u>Problem</u>	<u>Solution</u>
<b>The power indicator (refer to 4.1) is not lit.</b>	A. Check that the AC power cord is properly connected. B. Check the circuit breakers. If necessary, reset. C. Power input switch is in the off position
<b>The power indicator is lit, but Xenon lamp will not ignite.</b>	A. Check to make sure lamp cartridge is seated properly. B. Check that the lamp door is secured. C. Check hour meter; if lamp hours exceed the rated lamp life, replace the lamp (refer to 8.1).

## 11. CHART OF MEDICAL DEVICE SYMBOLS USED

	Manufacturer
	Date of manufacture (YYYY-MM-DD)
	"Authorized Representative" in the European Community
	Caution
	Caution: Federal law restricts this device to sale by or on the order of a licensed healthcare practitioner
	Consult Instructions for Use
	CE mark
	Not disposable in general waste,
	Caution: Hot Surface
	Caution: Dangerous Voltage
	Product Safety Mark
	Do not use if package is damaged

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