



LLS-2000 / LLS-1000 50W Solid State LED Lightsources



Premium LLS-2000 with 4-port turret



Standard LLS-1000 with 4-port turret

LLS-2000 / LLS-1000 Medical Grade 50W LED Lightsources

- Compact, efficient, durable sources of brilliant “daylight” quality light.
- Standard LLS-1000 provides brighter, whiter light than 150W Halogen units.
- Premium LLS-2000 utilizes new high efficiency design in same form, generating 40% more light (nominal) than LLS-1000 for more demanding applications.
- Color temperature: 6500° K nominal.
- Typical LED life: >50,000 hours.
- Low power consumption: 52W nominal.
- Electronic (PWM) intensity control from 0-100%.
- Four-port turret and single port versions available.
- Standard ports are ACMI, Olympus, Storz and Wolf; other ports available.
- Accepts 100-240V AC, 50/60 Hz input power for worldwide operation.
- Regulatory: UL60601-1, CAN/CSA C22.2 No. 601.1, EN60601-1-2; **CE** marked.
- Dimensions: 4.65” W x 3.32” H X 7.25” D (~10” overall).
- Weight : 2.5 lbs (1.14 kg).