

Headlights

CUDA[®]

SURGICAL

CUDA Surgical "RCS" Headlight

The Rear Cranial Support System introduces an innovative concept that greatly improves the comfort and experience of using a surgical headlight. Unlike the traditional side-to-side design, this system implements a crown support that extends from the front to the back, offering superior adjustability and stability for the module. By incorporating these advancements, user comfort is significantly enhanced, particularly during long and intricate procedures.



The Molded SureGrip Handle eases cable removal from the light source and helps prevent fiber damage.

CUDA Surgical Deluxe Headlight

Our Deluxe style headband is economical, lightweight, and suitable for daily use.

CUDA Surgical Sport Headlight

Our lightest headlight model, the SportBand provides superior comfort while ensuring module stability and adjustability



Surgical Headlight

Specifications

- A Lightweight Headband with Rear Cranial Support and Dual Ratchet adjustable Knobs for Premium Comfort
- Headlight Module with 20mm-110mm Variable Spot @ 16 inches/40cm working distance
- Fused 3.5mm Bifurcated Fiber Optic Cable 9' Long with Lensed Proximal End to optimize Light Output
- Bifurcated cable distributes weight evenly while molded SureGrip handle and reinforced sheathing protect the fiber bundle
- Headband with Module Weight - 4.5 oz/128 grams
- Headband, Module and Bifurcated Fiber Optic Cable Weight - 6.5 oz/180 grams

System Includes

- 1 Headlight With Rear Cranial Support "RCS"
- 1 20mm-110mm Variable Spot Module
- 1 3.5mm Bifurcated 9' Long SureGrip Cable with Lensed Proximal End
- 2 joysticks designed for precise positioning of the light spot - reusable and sterilizable
- 1 Duffel Bag for storage

Headlight Ordering Information

- DRH-110B-ESG: RCS Headband, 110mm variable spot module, 3.5 mm molded SureGrip bifurcated black cable
- DRH-110S-ESG: RCS Headband, 110mm variable spot module, 3.5 mm molded SureGrip single black cable
- SOH-110B-ESG: CUDA Elite SPORT Headband, 110mm variable spot module, 3.5 mm molded SureGrip bifurcated black cable
- SOH-110S-ESG: CUDA Elite SPORT Headband, 110mm variable spot module, 3.5 mm molded SureGrip single black cable
- ESG-3500B-BLACK: CUDA 3.5mm, molded SureGrip black bifurcated headlight cable, Wolf lensed
- ESG-3500S-BLACK: CUDA 3.5mm, molded SureGrip black single headlight cable, Wolf lensed
- SYC-0200: Rear cranial support headband
- SYC-0041: CUDA Surgical "RCS" only
- SYC0026: Module hinge extension for increased module mobility
- SYC-0040: Replacement comfort liners for headband (Set of 3)
- SYC-0014: Joysticks for use with all Cuda Headlights (5 per set, aluminum, autoclavable)

CUDA[®]
SURGICAL

6018 Bowdendale Avenue, Jacksonville, Florida 32216
USA Phone: 877-814-2237 International: +1 904-737-7611
www.CUDA.com

MADE IN THE USA

LIT 175B

CE



ISO 13485 Certified